

09 / 48 61 25

DATA ACQUISITION FOR MAGNETIC RESONANCE IMAGING TECHNIQUE

BACKGROUND OF THE INVENTION

FIELD OF THE INVENTION

The present invention relates to a method of data acquisition for a magnetic resonance imaging (MRI) technique, and more particularly, to a method of data acquisition for taking multiple images in sequence with a single contrast injection. According to the present invention, the multiple images can be taken with multi-orientation overlapping thin slabs, or the images can be individually taken as the patient is moved rapidly between different positions during the data acquisition period.

BACKGROUND AND SUMMARY OF THE INVENTION

Magnetic Resonance Imaging techniques for performing three-dimensional Magnetic Resonance Angiography (3D-MRA) are traditionally performed by placing the receiving antenna, or slab, over the relevant portion of the body, injecting the patient with a suitable dye and then triggering the M.R.I. data acquisition process while the patient holds his or her breath. Holding the breath is important to ensure that the arteries being examined do not move during data acquisition due to the expansion and contraction of the lungs and diaphragm. In conventional practice, the data are acquired during the arterial phase (before the dye is pumped by the heart into the venous system). Because the injected dye tends to dissipate quickly, in the arterial system, and is ultimately pumped into the venous system, data must be taken comparatively quickly. Thus, using conventional techniques, it has not been possible to image different regions of the body from a single dye injection. The time required to position the patient and imaging slab is too long and the injected dye will have dissipated before the next view can be set up and acquired.

When multiple regions must be imaged with a single dye injection, the conventional approach has been to use multiple coils positioned in advance over the regions of

interest. The data acquisition system is then programmed to collect data from the multiple coils in sequence. However, during the entire multiple slab acquisition sequence, the patient must remain very still and hold his or her breath. Unfortunately, for some tests, and for some critically ill patients, it is simply not possible to hold the breath throughout the entire sequence. (The data acquisition sequence for a single view can take approximately one to two minutes, depending on the data array size.)

The conventional data acquisition sequence involves the collecting of data for a single view and then processing that data to form an image suitable for display on screen or printout. A single view will typically comprise a plurality of partitions, with each partition comprising a plurality of lines. Typically, when imaging the chest, the partitions are acquired by electrocardiographic (E.C.G.) triggering at the R-wave. Thus, the number of partitions in a single view will correspond to the number of heartbeats. Within each heartbeat interval (R-R interval), a plurality of lines of data is acquired on the order of about ninety-six (96) lines. The number of lines acquired will depend on the patient's R-R interval and the time it takes to complete a single line of data acquisition.

In order to interpret the data meaningfully, certain parameters are set into the M.R.I. system. These parameters quantify such physical properties as the orientation of the coil, the pixel size, and the thickness of the slab or portion of the body to be imaged. The M.R.I. system uses these parameters to interpret the raw data from the receiving antenna coil, so that the data may be processed in a unified way by the imaging software. A typical view or slab might comprise twenty (20) to thirty (30) partitions, each being approximately two to three millimeters in thickness. In practice, a larger number of partitions cannot be completed in a reasonable, single

O
S
B
E
H
R
G
G
N
O
D
15
20
25

breath-hold period. However, with newer techniques, such as SYNC interpolation of data in the Fourier domain, the partitions can be doubled without increasing scan time.

After acquisition of the raw data, as described above, the M.R.I. system automatically converts this data into image data, using imaging software that is resident in the M.R.I. system. Image processing is very data intensive, accounting for a substantial portion of the overall data acquisition cycle. As noted above, the entire data acquisition cycle, including this processing of raw data, can take one to two minutes to complete. During this time, the injected dye will have dissipated or moved into the venous system, making subsequent views impractical from a single dye injection.

The present invention solves the aforementioned problems by altering the data acquisition cycle so that multiple views (in multiple orientations) can be taken. This is accomplished by delaying image processing until raw data from all views have been collected, and by programmatically inserting a predetermined dead time between views during which time the patient is permitted to exhale and inhale. The improved data acquisition system also automatically changes applicable parameters on the fly, to accommodate different orientations between views, or to alter other parameters as may be needed to obtain the most useful images. The ability to perform a sequence of views without the need for technician interaction between views is of considerable benefit in acquiring data from multiple slabs automatically. However, the improved data acquisition system does far more than this. In one embodiment, the data acquisition can supply position data to a moving table, allowing the patient to be automatically moved to different positions within the M.R.I. device. Accordingly, the patient can be moved using a kinematic examination table so that different body portions such as the chest, abdomen, and legs can be imaged in sequence with a single contrast injection.

094515
S12E5-O20
D20
094515
S12E5-O20
D20

Further areas of applicability of the present invention will become apparent from the detailed description provided hereinafter. It should be understood, however, that the detailed description and specific examples, while indicating preferred embodiments of the invention, are intended for purposes of illustration only, since various changes and modifications within the spirit and scope of the invention will become apparent to those skilled in the art from this detailed description.

10

BRIEF DESCRIPTION OF THE DRAWINGS

The present invention will become more fully understood from the detailed description and the accompanying drawings, wherein:

Figure 1 shows an example of the data acquisition stages, where three views are sequentially acquired at three different orientations;

Figure 2 illustrates diagrammatically how individual views are in turn subdivided into partitions, which are in turn subdivided into lines, a delay period is executed between each view;

Figures 3 and 4 give further details on how the views are subdivided into partitions (Fig. 4) and the partitions are subdivided into lines (Fig. 3);

25

Figure 5 shows how the software system of the improved data acquisition invention may be implemented; and

Figure 6 shows a data flow diagram for the control software according to the principles of the present invention.

DETAILED DESCRIPTION OF THE PREFERRED EMBODIMENTS

30

Figure 1 shows an example of the data acquisition system, where three views are sequentially acquired at three different orientations. Note that between each view is a dead time that is generated by the software of the invention. As illustrated, each view involves a number of physical parameter settings that the system uses in interpreting the raw data. As will be explained, the data acquisition system will automatically change parameters to

35

0015
0020
0025
0030
0035
0040
0045
0050
0055
0060
0065
0070
0075
0080
0085
0090
0095
0100
0105
0110
0115
0120
0125
0130
0135
0140
0145
0150
0155
0160
0165
0170
0175
0180
0185
0190
0195
0200
0205
0210
0215
0220
0225
0230
0235
0240
0245
0250
0255
0260
0265
0270
0275
0280
0285
0290
0295
0300
0305
0310
0315
0320
0325
0330
0335
0340
0345
0350
0355
0360
0365
0370
0375
0380
0385
0390
0395
0400
0405
0410
0415
0420
0425
0430
0435
0440
0445
0450
0455
0460
0465
0470
0475
0480
0485
0490
0495
0500
0505
0510
0515
0520
0525
0530
0535
0540
0545
0550
0555
0560
0565
0570
0575
0580
0585
0590
0595
0600
0605
0610
0615
0620
0625
0630
0635
0640
0645
0650
0655
0660
0665
0670
0675
0680
0685
0690
0695
0700
0705
0710
0715
0720
0725
0730
0735
0740
0745
0750
0755
0760
0765
0770
0775
0780
0785
0790
0795
0800
0805
0810
0815
0820
0825
0830
0835
0840
0845
0850
0855
0860
0865
0870
0875
0880
0885
0890
0895
0900
0905
0910
0915
0920
0925
0930
0935
0940
0945
0950
0955
0960
0965
0970
0975
0980
0985
0990
0995
1000
1005
1010
1015
1020
1025
1030
1035
1040
1045
1050
1055
1060
1065
1070
1075
1080
1085
1090
1095
1100
1105
1110
1115
1120
1125
1130
1135
1140
1145
1150
1155
1160
1165
1170
1175
1180
1185
1190
1195
1200
1205
1210
1215
1220
1225
1230
1235
1240
1245
1250
1255
1260
1265
1270
1275
1280
1285
1290
1295
1300
1305
1310
1315
1320
1325
1330
1335
1340
1345
1350
1355
1360
1365
1370
1375
1380
1385
1390
1395
1400
1405
1410
1415
1420
1425
1430
1435
1440
1445
1450
1455
1460
1465
1470
1475
1480
1485
1490
1495
1500
1505
1510
1515
1520
1525
1530
1535
1540
1545
1550
1555
1560
1565
1570
1575
1580
1585
1590
1595
1600
1605
1610
1615
1620
1625
1630
1635
1640
1645
1650
1655
1660
1665
1670
1675
1680
1685
1690
1695
1700
1705
1710
1715
1720
1725
1730
1735
1740
1745
1750
1755
1760
1765
1770
1775
1780
1785
1790
1795
1800
1805
1810
1815
1820
1825
1830
1835
1840
1845
1850
1855
1860
1865
1870
1875
1880
1885
1890
1895
1900
1905
1910
1915
1920
1925
1930
1935
1940
1945
1950
1955
1960
1965
1970
1975
1980
1985
1990
1995
2000
2005
2010
2015
2020
2025
2030
2035
2040
2045
2050
2055
2060
2065
2070
2075
2080
2085
2090
2095
2100
2105
2110
2115
2120
2125
2130
2135
2140
2145
2150
2155
2160
2165
2170
2175
2180
2185
2190
2195
2200
2205
2210
2215
2220
2225
2230
2235
2240
2245
2250
2255
2260
2265
2270
2275
2280
2285
2290
2295
2300
2305
2310
2315
2320
2325
2330
2335
2340
2345
2350
2355
2360
2365
2370
2375
2380
2385
2390
2395
2400
2405
2410
2415
2420
2425
2430
2435
2440
2445
2450
2455
2460
2465
2470
2475
2480
2485
2490
2495
2500
2505
2510
2515
2520
2525
2530
2535
2540
2545
2550
2555
2560
2565
2570
2575
2580
2585
2590
2595
2600
2605
2610
2615
2620
2625
2630
2635
2640
2645
2650
2655
2660
2665
2670
2675
2680
2685
2690
2695
2700
2705
2710
2715
2720
2725
2730
2735
2740
2745
2750
2755
2760
2765
2770
2775
2780
2785
2790
2795
2800
2805
2810
2815
2820
2825
2830
2835
2840
2845
2850
2855
2860
2865
2870
2875
2880
2885
2890
2895
2900
2905
2910
2915
2920
2925
2930
2935
2940
2945
2950
2955
2960
2965
2970
2975
2980
2985
2990
2995
3000
3005
3010
3015
3020
3025
3030
3035
3040
3045
3050
3055
3060
3065
3070
3075
3080
3085
3090
3095
3100
3105
3110
3115
3120
3125
3130
3135
3140
3145
3150
3155
3160
3165
3170
3175
3180
3185
3190
3195
3200
3205
3210
3215
3220
3225
3230
3235
3240
3245
3250
3255
3260
3265
3270
3275
3280
3285
3290
3295
3300
3305
3310
3315
3320
3325
3330
3335
3340
3345
3350
3355
3360
3365
3370
3375
3380
3385
3390
3395
3400
3405
3410
3415
3420
3425
3430
3435
3440
3445
3450
3455
3460
3465
3470
3475
3480
3485
3490
3495
3500
3505
3510
3515
3520
3525
3530
3535
3540
3545
3550
3555
3560
3565
3570
3575
3580
3585
3590
3595
3600
3605
3610
3615
3620
3625
3630
3635
3640
3645
3650
3655
3660
3665
3670
3675
3680
3685
3690
3695
3700
3705
3710
3715
3720
3725
3730
3735
3740
3745
3750
3755
3760
3765
3770
3775
3780
3785
3790
3795
3800
3805
3810
3815
3820
3825
3830
3835
3840
3845
3850
3855
3860
3865
3870
3875
3880
3885
3890
3895
3900
3905
3910
3915
3920
3925
3930
3935
3940
3945
3950
3955
3960
3965
3970
3975
3980
3985
3990
3995
4000
4005
4010
4015
4020
4025
4030
4035
4040
4045
4050
4055
4060
4065
4070
4075
4080
4085
4090
4095
4100
4105
4110
4115
4120
4125
4130
4135
4140
4145
4150
4155
4160
4165
4170
4175
4180
4185
4190
4195
4200
4205
4210
4215
4220
4225
4230
4235
4240
4245
4250
4255
4260
4265
4270
4275
4280
4285
4290
4295
4300
4305
4310
4315
4320
4325
4330
4335
4340
4345
4350
4355
4360
4365
4370
4375
4380
4385
4390
4395
4400
4405
4410
4415
4420
4425
4430
4435
4440
4445
4450
4455
4460
4465
4470
4475
4480
4485
4490
4495
4500
4505
4510
4515
4520
4525
4530
4535
4540
4545
4550
4555
4560
4565
4570
4575
4580
4585
4590
4595
4600
4605
4610
4615
4620
4625
4630
4635
4640
4645
4650
4655
4660
4665
4670
4675
4680
4685
4690
4695
4700
4705
4710
4715
4720
4725
4730
4735
4740
4745
4750
4755
4760
4765
4770
4775
4780
4785
4790
4795
4800
4805
4810
4815
4820
4825
4830
4835
4840
4845
4850
4855
4860
4865
4870
4875
4880
4885
4890
4895
4900
4905
4910
4915
4920
4925
4930
4935
4940
4945
4950
4955
4960
4965
4970
4975
4980
4985
4990
4995
5000
5005
5010
5015
5020
5025
5030
5035
5040
5045
5050
5055
5060
5065
5070
5075
5080
5085
5090
5095
5100
5105
5110
5115
5120
5125
5130
5135
5140
5145
5150
5155
5160
5165
5170
5175
5180
5185
5190
5195
5200
5205
5210
5215
5220
5225
5230
5235
5240
5245
5250
5255
5260
5265
5270
5275
5280
5285
5290
5295
5300
5305
5310
5315
5320
5325
5330
5335
5340
5345
5350
5355
5360
5365
5370
5375
5380
5385
5390
5395
5400
5405
5410
5415
5420
5425
5430
5435
5440
5445
5450
5455
5460
5465
5470
5475
5480
5485
5490
5495
5500
5505
5510
5515
5520
5525
5530
5535
5540
5545
5550
5555
5560
5565
5570
5575
5580
5585
5590
5595
5600
5605
5610
5615
5620
5625
5630
5635
5640
5645
5650
5655
5660
5665
5670
5675
5680
5685
5690
5695
5700
5705
5710
5715
5720
5725
5730
5735
5740
5745
5750
5755
5760
5765
5770
5775
5780
5785
5790
5795
5800
5805
5810
5815
5820
5825
5830
5835
5840
5845
5850
5855
5860
5865
5870
5875
5880
5885
5890
5895
5900
5905
5910
5915
5920
5925
5930
5935
5940
5945
5950
5955
5960
5965
5970
5975
5980
5985
5990
5995
6000
6005
6010
6015
6020
6025
6030
6035
6040
6045
6050
6055
6060
6065
6070
6075
6080
6085
6090
6095
6100
6105
6110
6115
6120
6125
6130
6135
6140
6145
6150
6155
6160
6165
6170
6175
6180
6185
6190
6195
6200
6205
6210
6215
6220
6225
6230
6235
6240
6245
6250
6255
6260
6265
6270
6275
6280
6285
6290
6295
6300
6305
6310
6315
6320
6325
6330
6335
6340
6345
6350
6355
6360
6365
6370
6375
6380
6385
6390
6395
6400
6405
6410
6415
6420
6425
6430
6435
6440
6445
6450
6455
6460
6465
6470
6475
6480
6485
6490
6495
6500
6505
6510
6515
6520
6525
6530
6535
6540
6545
6550
6555
6560
6565
6570
6575
6580
6585
6590
6595
6600
6605
6610
6615
6620
6625
6630
6635
6640
6645
6650
6655
6660
6665
6670
6675
6680
6685
6690
6695
6700
6705
6710
6715
6720
6725
6730
6735
6740
6745
6750
6755
6760
6765
6770
6775
6780
6785
6790
6795
6800
6805
6810
6815
6820
6825
6830
6835
6840
6845
6850
6855
6860
6865
6870
6875
6880
6885
6890
6895
6900
6905
6910
6915
6920
6925
6930
6935
6940
6945
6950
6955
6960
6965
6970
6975
6980
6985
6990
6995
7000
7005
7010
7015
7020
7025
7030
7035
7040
7045
7050
7055
7060
7065
7070
7075
7080
7085
7090
7095
7100
7105
7110
7115
7120
7125
7130
7135
7140
7145
7150
7155
7160
7165
7170
7175
7180
7185
7190
7195
7200
7205
7210
7215
7220
7225
7230
7235
7240
7245
7250
7255
7260
7265
7270
7275
7280
7285
7290
7295
7300
7305
7310
7315
7320
7325
7330
7335
7340
7345
7350
7355
7360
7365
7370
7375
7380
7385
7390
7395
7400
7405
7410
7415
7420
7425
7430
7435
7440
7445
7450
7455
7460
7465
7470
7475
7480
7485
7490
7495
7500
7505
7510
7515
7520
7525
7530
7535
7540
7545
7550
7555
7560
7565
7570
7575
7580
7585
7590
7595
7600
7605
7610
7615
7620
7625
7630
7635
7640
7645
7650
7655
7660
7665
7670
7675
7680
7685
7690
7695
7700
7705
7710
7715
7720
7725
7730
7735
7740
7745
7750
7755
7760
7765
7770
7775
7780
7785
7790
7795
7800
7805
7810
7815
7820
7825
7830
7835
7840
7845
7850
7855
7860
7865
7870
7875
7880
7885
7890
7895
7900
7905
7910
7915
7920
7925
7930
7935
7940
7945
7950
7955
7960
7965
7970
7975
7980
7985
7990
7995
8000
8005
8010
8015
8020
8025
8030
8035
8040
8045
8050
8055
8060
8065
8070
8075
8080
8085
8090
8095
8100
8105
8110
8115
8120
8125
8130
8135
8140
8145
8150
8155
8160
8165
8170
8175
8180
8185
8190
8195
8200
8205
8210
8215
8220
8225
8230
8235
8240
8245
8250
8255
8260
8265
8270
8275
8280
8285
8290
8295
8300
8305
8310
8315
8320
8325
8330
8335
8340
8345
8350
8355
8360
8365
8370
8375
8380
8385
8390
8395
8400
8405
8410
8415
8420
8425
8430
8435
8440
8445
8450
8455
8460
8465
8470
8475
8480
8485
8490
8495
8500
8505
8510
8515
8520
8525
8530
8535
8540
8545
8550
8555
8560
8565
8570
8575
8580
8585
8590
8595
8600
8605
8610
8615
8620
8625
8630
8635
8640
8645
8650
8655
8660
8665
8670
8675
8680
8685
8690
8695
8700
8705
8710
8715
8720
8725
8730
8735<br

match the physical attributes of each given view or orientation. The different orientations can include taking an image of, for example, a vertical slab or slice of a patient's chest, and then taking images at a different orientation such as a horizontal slab or slice of the same area giving an orientational view which is different than the first image. Most importantly, according to the present invention, a delay period is provided between the first imaging sequence and the second image sequence in order to allow the patient to exhale and inhale prior to again holding his breath for the second image to be taken. In addition, the patient can be moved on a kinematic imaging table during the delay period so that different portions of the patient's body can be imaged in sequence using a single injection of contrast material as the contrast material travels through the vascular system to the different portions of the patient's body. The system of the present invention thus permits a dynamic contrast-enhanced breath-hold magnetic resonance angiography technique using a multi-orientation and multiple overlapping slab acquisition (MO-MOTSA), as well as an MRA technique which allows multiple body portions to be scanned in sequence using a single contrast injection.

Figure 2 illustrates diagrammatically how individual views are in turn subdivided into partitions, which are in turn subdivided into lines. Note that after all partitions of a given view have been acquired, the system will automatically set a predetermined delay before the partitions for the next view are acquired. Figures 3 and 4 give further details on how the views are subdivided into partitions (Fig. 4) and partitions are subdivided into lines (Fig. 3). The manner in which views are subdivided into partitions and lines are generally known in the art.

Figure 5 shows how the software system of the improved data acquisition invention may be implemented. Although there are a number of different ways to achieve the data acquisition sequence, the illustrated embodiment uses a

5
1000435225.01200
2025
30

35

conditional branching loop to cycle through the plurality of views, selecting parameters from a parameter table. In Figure 5, the parameter table is shown diagrammatically at 100. In the illustrated example, there are four parameters (a, b, c, d) for each of three views (1, 2, 3). The parameters may be supplied by user input through a suitable interface 102. If desired, multiple parameter tables can be stored in the system and used as templates for performing different types of M.R.I. techniques. The user would then select the appropriate template corresponding to the type of procedure being performed, and would then enter into the template the desired parameters corresponding to the patient's physical characteristics or to the radiologist's instructions.

The software routine begins at Step 110 by pointing to the head of the parameter table. As will be described, the pointer is used to determine what parameters are loaded from the table. The pointer is indexed to the next set of parameters in the table for each of the views to be performed.

After initializing to the head of the parameter table, the routine, in Step 112, loads the parameters from the table into the system for execution. Then, in Step 114, after the patient has been injected with a suitable contrast material and the material has traveled through the arterial system to the body portion being viewed, the scan for the current view is run. This will involve acquiring all lines of data for all partitions that comprise the current view. See Figures 2-4 for details. After the scan for the current view is complete, the system then executes a dead time routine at 116. Preferably, the dead time interval can be a predetermined value supplied by the user through user interface 102. Although the current system uses the same dead time value between all views, the system can be readily adapted to execute different dead times between views, as may be desired in a particular procedure. Typically, the dead time will correspond to the amount of

D0115
S125
E25
-020
D000

time it takes for a patient to exhale and inhale and again hold his or her breath for the next view to be taken at a different orientation. When the patient is being moved so that a different body portion (such as moving from the abdomen to the legs) is being imaged, the dead time can correspond to the time it takes for the contrast material to travel from the first body portion being viewed to a second body portion to be viewed.

It will be recalled that during the dead time the patient is permitted to exhale and inhale. If desired, the system can be programmed to issue a prompt at Step 118 to notify the attending technician or patient that the patient may exhale and inhale. A suitable audible or visual prompt can be generated for this purpose. Also, if the system is configured to perform automatic patient repositioning, such as on a kinematic table, the appropriate motor control signals are generated to move the table on which the patient is resting at Step 120. The step of moving the patient can also be performed manually by a technician.

After the dead time is complete, the program tests at conditional branching operation 122 whether the final scan specified in the parameter table has been completed. If not, control branches to Step 124, where the parameter table pointer is indexed to the next entry in the table, and the program branches back to Step 112. At this time, additional images are obtained in sequence until the desired orientations and/or the desired portions of the body are imaged in series. If the final scan has been completed at Step 122, then control branches to Step 126, where the image data for all views are then processed.

Figure 6 illustrates a data flow diagram for the control system software according to the principles of the present invention. Existing MRI systems include a data collect module 200 which loads the parameters for a scan sequence and collects the image data during the scan process. The MRI systems also include an imaging processing module 202 which receives the image data from

the data collect module 200 and processes the images and displays them either on computer screen or in hard copy form as represented by reference numeral 204. Existing MRI systems also include means for entering and storing the parameters 206 for a scan sequence, as well as storage means 208 for storing the image data collected by the data collect module 200. The image processing module 202 receives the image data from the data storage 208. Typically, MRI systems are utilized for a single scan sequence and the data collected by the data collect module 200 is immediately processed by the image processing module 202 upon completion of the scan sequence. Existing MRI systems are not provided with software for carrying out a plurality of scan sequences using different sets of parameters which are separated by a delay period and storing the image data corresponding to each scan sequence in designated locations so that the image processing module can, upon completion of the scan sequences, process the image data for each of the scan sequences (views).

Accordingly, the present invention provides a controller 210 which initiates each scan sequence and has associated with it a parameters table manager 212 which manages the parameter sets, stored in the parameters table 206 which are input into the data collect module 200 corresponding to each scan sequence. The controller also has associated with it a data store manager 214 for managing the storage of image data corresponding to each scan sequence. The data store manager 214 tells the data collect module 200 where to store the image data which is collected by the data collect module 200 and stored in storage 208.

Controller 210 also has associated with it, a delay timer 216 which institutes a delay period between each scan sequence. The controller 210 initiates a first scan sequence by directing the parameters table manager 212 to input the appropriate parameters from parameter tables 206 into the data collect module 200. The controller also

directs the data storage manager 214 to instruct the data collect module 200 where to store the image data collected in the scan sequence. At the completion of a scan sequence, the controller initiates the delay timer 216 for instituting a delay period between each scan sequence. During the delay period, the patient is allowed to breathe prior to holding his or her breath again during a subsequent scan sequence at a different orientation. The patient can also be moved to a second location so that a second scan sequence can be performed on a second body portion. If an automatic table moving device is utilized, such as table motor 218, controller 210 can provide a signal to actuate the table motor 218 to move the patient to a second orientation for obtaining a scan sequence of a second body portion of the patient. A control interface 220 is provided for activating the image processing module 202 in response to a signal from the controller 210 after the data acquisition is complete. The software system of the present invention allows a method of data acquisition for taking multiple images in sequence with a single contrast injection wherein the images can be at different orientations (MO-MOTSA) or wherein the patient is moved rapidly between different positions during the data acquisition. In either case, a delay period is provided between each scan sequence to allow the patient to breathe and again hold his or her breath for a subsequent scan sequence.

Description of Clinical Examples

Three dimensional magnetic resonance angiography (3D-MRA) techniques have been routinely used in the past for vascular imaging. In particular, 3D-MRA in a single breath-hold [Simonetti OP, Finn JP, White RD, Bis KG, Shetty AN, et al. ECG-triggered breath-held gadolinium-enhanced 3D MRA of the thoracic vasculature(abstr). In: Book of abstract: International society of Magnetic Resonance in Medicine 1996:703. Shetty AN, Shirkhoda A, Bis KG, Alcantara A. Contrast enhanced three dimensional MR

angiography in a single breath hold: A novel technique. AJR 1995;165:1290-1292.] has been shown to be a very useful and rapid technique when used with paramagnetic contrast agents. Given the rapid nature of this new sequence, with the acquisition of an entire 3D-data set during the passage of MR contrast material through the arterial system, arterial phase imaging is possible. A delayed acquisition following the initial arterial phase results in the visualization of venous flow. Conventionally, a breath-hold 3D-MRA study is performed by placing a single slab (view) of a desired thickness whose orientation is set to cover most of the vessels in question. The number of partitions within the slab (view) and in-plane phase encoding steps are adjusted to complete the acquisition in a reasonable single breath-hold period which is anywhere from 15-25 seconds. The data processing takes anywhere from 20-40 seconds depending on the data array size. These techniques have been successfully used to obtain a dynamic measurement in a single orientation. [Maki JH, Prince MR, Londy F, Chenevert TL. The effect of time varying intravascular signal intensity and k-space acquisition order on 3D MR angiography image quality (abstr). In: Book of abstract: International Society of Magnetic Resonance in Medicine 1996:237. Debatin JF, Schmidt M., Gohde S., Hany T., Pfamatter T, Krestin GP., McKinnon GC. 3D breath hold MRA of the renal arteries under 30 seconds. (abstr). In: Book of abstract: International Society of Magnetic Resonance in Medicine 1996:165.] In many cases, multiple measurements are made in the same orientation to obtain arterial and venous phase information. However, to obtain additional images at a different orientation, the data acquisition is repeated by repositioning the slab (view) with a new orientation. A different orientation is frequently required to better portray the 3-dimensional aspects of complex vascular anatomy with high resolution. Since, the data acquisition of the newly positioned slab (view) can only commence after the previous image data

DRAFT
15
20
25
30
35

processing is completed, a wait of about 1-2 minutes between scans is elapsed during which most of the contrast bolus is washed out. A subsequent measurement will yield images with poor arterial and venous signal due to this prolonged delay. The present invention provides a method in which multiples slabs, either in the same or different orientation may be prescribed to acquire the arterial and venous phases with a single contrast injection.

Materials and methods:

All studies were performed using a 1.5 Tesla Magnetom Vision MR system (Siemens Medical Systems, Iselin, NJ). Thirty subjects with known aortic (10 thoracic, 1 abdominal, 1 arch vessel), pulmonary artery (6), or cardiac (7) disease and 5 volunteers were recruited with informed consent approved by the institutional review board of the hospital. Patients were positioned supine with a quadrature phased array body coil. A 20-gauge angio-cath needle was placed in an antecubittal vein and the other end of the line was connected to a contrast injector (Spectra II, MedRad, Pittsburgh, PA) which was used to infuse contrast at a rate of 2 or 3 cc /sec. Any 3D MRA or MRI pulse sequence can be used with this technique. We used an ECG triggered 3D FLASH (Fast Low Angle single Shot) pulse sequence. Each of the partitions along the slab (view) select direction was acquired by triggering at the R-wave. Thus, the number of partitions correspond to the number of heart beats. About 96 lines were acquired in a single R-R interval for each partition. The number of lines that were acquired depended on the patient's R-R interval and the time to complete a single line of data acquisition. The imaging parameters were TR/TE = 5.0/2.0, flip angle = 15°, matrix = 96 x 256, partition thickness = 2 - 3 mm and the number of partitions were limited to 22-30 so that a single measurement could be completed in a reasonable single breath hold (after deep inspiration) period. It should be noted that the above sequence is just one sequence that was implemented. Other gradient echo, spin echo, or fast

(turbo) spin echo sequences can also be implemented using this concept.

The delay between successive slab measurements was fixed at eight seconds. During this period, patients were prepared to hold their breath for the next measurement. Prior to contrast injection, the sequence was applied once to obtain base-line images in the same three orientations which were then used for subtracting from the corresponding post-contrast images. Subtraction was used in only five patients whereas in remaining patients source images resulting from contrast were used. Prohance (Gadoteridol) (Bracco Diagnostics, Princeton, NJ) at the dose of 0.2 mmol/kg was injected at the rate of 2 or 3 cc/sec using the injector followed by 15 cc of saline flush. The first measurement was obtained between 5-10 seconds (for pulmonary studies) or 10-15 seconds (for aortic studies) after the onset of contrast injection. Second and third measurements were subsequently obtained with a total time for all three measurements of 54 seconds.

Results:

All patients tolerated three successive measurements with breath-holding. All measurements revealed a moderate enhancement of the blood pool of the target area of interest. All 30 subjects underwent three measurements in which 26 were imaged in two different planes (one of the imaging planes was repeated) over time and four were imaged with three different planes over time. Among the 10 patients with thoracic aortic disease, the following abnormalities were imaged: Type A dissections (4), Type B dissections (2), aorto-annular ectasia (1), atherosclerotic aneurysms (2) and atherosclerotic narrowing (1). There were two normal thoracic aortas imaged. One plane of imaging is not optimum in dissection patients since the imaging plane may be parallel with the intimal flap which may go undetected. The three-dimensional relationships of thoracic aortic aneurysms was also better appreciated with the MO-MOTSA technique.

09435
20
21
22
23
24
25
26
27
28
29
30
31
32
33
34
35

Among the pulmonary abnormalities, the following were imaged with MO-MOTSA: pulmonary artery aneurysm (1), lung carcinoma mimicking arteriovenous malformation on CT (1), infiltrate/atelectasis (2) and chronic pulmonary embolus (1).

5

There were seven patients with a variety of cardiac abnormalities that were imaged with MO-MOTSA and included the following: hydatid cyst of the interventricular septum (1), pericardial effusions (3), pericardial cyst (1), right atrial lipoma (1) and left ventricular aneurysm with thrombus (1). Thickening of the pericardium with pericardial effusions was better appreciated with the multi-planar aspects of MO-MOTSA imaging.

DRAFTED
REVIEWED
APPROVED
DISTRIBUTED
20

Discussion:

We have shown a new way of acquiring 3D-MRA images in multiple orientations. The technique can also be used in other parts of the body for depicting vascular and non-vascular anatomy. Typically, a 3D pulse sequence takes 15-18 seconds to complete 96 lines and 24 partitions. By running the pulse sequence repeatedly, arterial and delayed arterial and venous phase imaging is possible. It should be noted that other sequences can also be implemented. Most scanners provide a means to achieve this by increasing the number of measurements of the same slab orientation with an appropriate time delay between measurements. However, this restricts one to achieving vascular information along only one orientation. To render projections along the other directions, a multi-planar reformatting (MPR) of the data can be used before subjecting the data to maximum intensity projection (MIP), however, the spatial resolution is lost especially when the raw data voxels are non-isotropic. In addition, with MIP, due to pixel replication, the projections deviating significantly from the true acquisition plane have a significantly diminished resolution. The method described here allows for acquiring data in multiple orientations without a loss of vascular signal. This also allows for

25

30

35

MIP along those multi-planar directions without the loss of spatial resolution.

Using a variable delay between 3D slabs allowed for imaging of vascular territories in different orientations during the arterial peak of contrast enhancement. The method is very useful in obtaining high resolution images of the aorta and pulmonary vasculature and cardiac anatomy. This method may also prove useful in the future when the patient table is moved to center different anatomic locations in the center of the magnet when performing aortic runoff evaluations or when the need for imaging the thoracic and abdominal aorta arises.

The invention being thus described, it will be obvious that the same may be varied in many ways. Such variations are not to be regarded as a departure from the spirit and scope of the invention, and all such modifications as would be obvious to one skilled in the art are intended to be included within the scope of the following claims.